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User Guide

www.wtavet.com.br

Incubadora de Bancada EVE

EVE Benchtop Incubator

Eve[®]

ADVERTÊNCIA

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Index

WARNING.....	26
Copyright.....	26
WELCOME.....	27
GENERAL DESCRIPTION.....	27
ACCESSORIES.....	28
TECHNICAL SPECIFICATIONS	28
IMPORTANT PRECAUTIONS	29
Repairs.....	29
Transportation	29
Cleaning.....	29
Do Not Set Items on Incubator	29
Power Cord.....	29
Electrical Surges.....	29
INPUTS AND POWER BUTTON.....	30
POWER SUPPLY	30
ASSEMBLY OF THE GAS PREHEATING SYSTEM.....	30
CONNECTING THE GAS PNEUMATIC HARNESS	31
INCUBATOR START UP	31
INITIAL SET UP OF THE INCUBATOR.....	33
Gas Inlet Regulation.	33
Adjusting the Chamber Temperature.....	33
Adjusting the Chamber Pressure	34
Adjusting the Time and Date	34
PREPARING CELL CULTURE CHAMBERS	35
AUTOMATIC REPLACEMENT OF GAS ENVIRONMENT	35
HANDLING DURING CULTURE.....	36
CLEANING CULTURE CHAMBER.....	36
ALERTS	37
Wi-Fi	37
REMOTE CONTROL (WTA WEBSITE)	38
Push Notifications	43
WARRANTY	44
PROBLEM RESOLUTION	46
TECHNICAL ASSISTANCE / CUSTOMER SERVICE	47

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Figure 4 consists of two screenshots labeled A and B. Screenshot A shows the 'Registrar-se' (Register) form. It includes fields for 'Nome de Usuário' (User Name), 'E-mail' (Email), 'Senha' (Password), 'Confirmar Senha' (Confirm Password), 'Chave de Acesso' (Access Key) with the note 'Fornecida na compra' (Provided at purchase), and a checkbox for accepting 'Termos de Privacidade' (Privacy Terms). Screenshot B shows the 'Conecte-se' (Connect) form. It includes fields for 'E-mail' and 'Senha' (Password), and a link for 'Esqueceu a senha?' (Forgot password?). Both forms have a red 'Next Step' button at the bottom right.

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The screenshot shows a mobile application interface for monitoring three chambers. At the top, it displays "EVE18014" and a bell icon. Below is a table with the following data:

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At the bottom are three buttons: "CONFIGURAÇÕES" (with a gear icon), "DATALOG" (with a document icon), and "REMOVER" (with a delete icon).

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e-mail: sac@wtavet.com.br

Warning

To avoid accidents and/or damage, carefully read and follow the instructions contained in this manual.

Symbol indicates
electrical shock hazard

Symbol indicates important
operating information

Copyright

EVE Benchtop Incubator User's Manual.

WTA, Cravinhos, São Paulo, Brazil.

The contents of this manual are protected by copyright laws and may not be reproduced, distributed, disseminated or stored without the prior written permission of WTA - Watanabe Applied Technology - EPP.

For authorization, additional information and questions, please contact our technical assistance department by email, letter or telephone.

Service from Monday to Friday, 8:30 am to 5:30 pm (GMT -03:00)

Telephone: + 1 (979) 324-6168

e-mail: usa1@wtatech.com

Comments

The information, illustrations and specifications included in this manual describe the specific features of the current EVE benchtop incubator model.

Due to development of new technologies and raw materials, WTA reserves the right to change the design and technical specifications of the product whenever necessary and without prior notice, thus not incurring obligations of any kind.

This manual is an integral part of the product. Therefore, the consumer should keep it for reference.

WELCOME

Congratulations on your decision to purchase a WTA product, a brand that guarantees quality and technology in all products that it manufactures.

WTA is a Brazilian company specializing in solutions for advanced animal reproduction technologies. Our products are developed to work with cattle, sheep and horses, recognized by the most renowned experts, veterinarians and scientists in the world.

Its anodized aluminum construction offers, in addition to an elegant design, internal and external surfaces that are easy to clean, durable and corrosion resistant.

Please read manual thoroughly before using the product.

GENERAL DESCRIPTION

Essential for cell culture, incubators reproduce the ideal physiological conditions for cell growth and development.

The EVE benchtop incubator has an electronic control system that keeps conditions, including temperature and gas pressure, extremely precise. Its small -size cultivation chambers, combined with an efficient sealing system, provide quick recovery and maintenance of the gas conditions after opening and closing its lids. These features allow you to safely carry out oocyte maturation, in vitro fertilization, embryo and general cell culture.

Its anodized aluminum and stainless steel construction offers a durable and decorative finish that is easy to clean to maintain the most hygenic surfaces for culture conditions.

This incubator has microprocessors and a 4.3" color liquid crystal screen with 480 x 272 pixels (TFT-LCD "Thin Film Transistor liquid crystal display") that allow for user friendly touch screen programming.

The USB "bootloader" function makes it possible to easily access and update WTA incubator programs.

Note: Incubator models available with 1 (modular) or 3 incubation chambers

Note: The incubator can be customized with a "high lid" in some or all cultivation chambers. The high lid provides a 60% increase in internal volume and accommodates WTA culture microtubes vertically.

ACCESSORIES

When unpacking your equipment, make sure your accessories are included in the package:

- 1 Power cord.
- 1 Pneumatic harness with quick connect valve connectors.
- 1 Electrical power supply.
- 3 (extra) pairs of pre-heating system hoses *.
- 3 Pairs (extra) of sealing “rings” for cork stoppers *.
- 8 Aluminum racks for micro culture tubes.
- 3 Gas preheating systems *.

* Modular model: only one extra pair and one pre-heating system.

Note: If any accessory is missing or damaged, please contact the WTA customer service department within 3 days from the date of delivery of the product.

TECHNICAL SPECIFICATIONS

- Dimensions (LxWxH): 22.2x15.7x7.1 in. Modular: 8.3x15.7x7.1 in.
- Weight: 22 kg. Modular: 8 kg
- Maximum power: 200W. Modular: 90 W.
- Equipment DC supply: 12 V – 16.66A. Modular: 12Vdc – 7.5A.
- Power AC Supply: 100-240Vac, 3~1A, 50/60hz.
- Temperature range: 34 - 45° C (93-113 F).
- Chamber pressure range: 1 - 50 mmHg.
- Visual alarms on display screen
- Chamber temperature accuracy: ± 0.5°C (1 F).
- Gas inlet pressure range: 9 to 11 Psi (0.6kgf/cm² or 0.8 kgf/cm²).

IMPORTANT PRECAUTIONS

Repairs:

Do not open the equipment casing. All repairs and maintenance must be carried out by WTA technical assistance service or authorized professionals.

Transportation:

Avoid vibrations during transport which could damage the incubator. Care should be taken to properly pad and protect the equipment against accidents.

Cleaning:

Clean the equipment with a soft cloth slightly dampened with water. Do not use alcohol or other solvents as they will damage the finish of the equipment.

Do Not Set Items on Incubator:

Do not place liquid or coins, paper clips or other metallic objects on the control module. If liquid or any objects drops inside the equipment, turn it off immediately and contact technical assistance.

Power Cord:

Do not step on, apply weight to, or kink the power cord.

Electrical Surges:

It is recommended to use a surge protector. In the event of thunderstorms with electrical discharges, turn off the incubator and unplug .

INPUTS AND POWER BUTTON

The rear panel (Figure 1) includes the power buttons (1), USB inputs (2), power supply input (3), gas mixture input (4) and gas chambers output (5).


Fig. 1 - Rear panel details.

Note. The power buttons (1) located on the rear panel is used to completely stop the operation of each individual chamber. They must be kept in the off position when the equipment is not being used.

POWER SUPPLY

The power supply cord powers all three chambers of the incubator and provides some protection against power fluctuations and surges that may occur in the electrical network.

ASSEMBLY OF THE GAS PREHEATING SYSTEM

Connect the end of the smaller diameter hose that has the filter (1) to the aluminum stopper attached to the conical tube. Then connect the other end to the port located at the base of the Incubator chamber (Figure 2). Connect one end of the larger diameter hose (2) to the other connector on the conical tube stopper and the other end of the hose to the port located on the Incubator lid. Add 20 mL of filtered water to the conical tube ("Bubbler"),

tightly screw on the aluminum stopper and insert the assembly into the tube holder. When necessary, clean the inside of the conical tube with a solution of water and a neutral detergent.


Fig. 2 - Gas preheating system set up.

Note: Regularly check the water level in the tube and change/replace water weekly.

CONNECTING THE GAS PNEUMATIC HARNESS

The pneumatic hose connectors must be coupled to the gas mixture inlets of the chambers (Fig. 1 - 4) and the other end to a pressure regulator (item not supplied). WTA recommends using a pressure regulator (manometer) with the following characteristics: Maximum inlet pressure 220kgf/cm²; Maximum outlet pressure = 3kgf/cm² and Maximum flow = 30m³/h).

The regulator must be attached to the gas cylinder containing the gas mixture.


Keep the pressure regulator CLOSED until the incubator is turned on, then the inlet pressure adjustment can be made.

The connectors can be disconnected from the incubator without leaking or compromising the other chambers as they have a gas lock system.

INCUBATOR START UP

After installing the hoses, connect the power supply cable to the rear connector of the incubator(Fig.1- 3). Press the power button (Fig.1- 1) to the ON position. The red “LED” light

inside the button will illuminate and the display will show the initial “Home” screen, and the touchscreen will become functional.


MENU: Access to the following options:


Home: Return to the home screen

Operation: Access operations related to: Gas, WiFi, Info, and USB input.

Config: Access the cultivation parameters configuration screen.

Service: Restricted to use by WTA service

The date, time, temperature, and extra pressure settings (Config) must be assigned during initial set up before first use of the equipment.

INITIAL SET UP OF THE INCUBATOR

Gas Inlet Regulation

After turning on the Incubator, the home screen will be displayed and will show "Inlet Pressure" on the right side of the display. Open the cylinder outlet pressure regulator valve until the pressure gauge reads 0.70 kgf/cm². Observe "Inlet Pressure" on the incubator display and make small adjustments (by turning the regulator) until you get an inlet pressure level equal to 10 PSI.

Outlet pressure greater than 0.9 kgf/cm² can damage equipment components which are not covered under warranty.

If the pneumatic hose connector is not connected to the incubator, the initial screen will show the following image which will flash.


In the configuration screen ("Config") the following settings can be adjusted:

- date and time
- chamber pressure (1 to 50 mmHg),
- incubation chamber temperature (34 to 45⁰ C/ 93 to 113⁰ F),
- temperature of the conical tube of the gas mixture pre-heating system ("Bubbler").

Adjusting the Temperature


Select "Menu" => "Config" => "Temperature".

Use "+" or "-" to adjust to the desired temperature, then press => "Confirm".

Adjusting the Chamber Pressure

Select "Menu" => "Config" => "Pressure".


Use "+" or "-" to adjust to the desired pressure, then press => "Confirm".

NOTE: The minimum value should be 1 mmHg.

Adjusting the Time and Date

Select "Menu" => "Config" => "Date/Clock"

Click in the box indicated and enter the values for Year- Month- Day: Hour- Minutes- Seconds consecutively, then press "Ok". "Esc" will cancel the date/ time input.

Once the temperature, pressure, date and time values are adjusted, the equipment will initiate temperature and pressure adjustments of the incubation chambers.

Note: Wait until the internal temperature has stabilized before adding culture plates or tubes.

Display: Adjusts the brightness level of the display

Adjust to the desired brightness and hold for 2 seconds. Press "back" to return to the previous screen .

PREPARING CELL CULTURE CHAMBERS

After the initial set up, cell culture dishes can be introduced into the incubation chamber. On the bottom plate of the culture chamber (Fig. 3) there are demarcations to help arrange the plates. The humidity reservoir (1) should be filled with filtered water while the incubator is in use and be placed on the central gas exhaust pin. Do not allow reservoir to dry while during incubations.

Avoid overflowing to prevent water from entering the gas exhaust holes. Liquid could damage the chamber control valves.


Fig. 3 - Bottom plate of the culture chamber


Upon closing the lid of the cultivation chamber, the gas environment will be adjusted automatically. This process will occur every time the cover is opened and closed again to maintain culture equilibrium.

If a gas leak is detected, the gas exchange process will be stopped and a visual alarm message will be displayed on the screen.


AUTOMATIC RECIRCULATION OF GAS ENVIRONMENT

The exchange intervals (minutes) of the gas conditions and the gas circulation time (seconds) inside the cultivation chambers can be programmed depending on culture requirements.

Select: "Menu" =>"Operation" =>"Gas".

Select "Yes" and then use "+" or "-" to adjust the timing.

Then press => "Confirm".

The equipment will replace the high O₂ atmosphere by releasing small volumes of the gas mixture until reaching the set pressure and temperature.

When automatic gas change is enabled the "home" screen will display the following image: **GAS CIRCULATION** . If the option is not enabled, the screen will display: **GAS CIRCULATION** .

ACTIVATED

DEACTIVATED

HANDLING DURING CULTURE

The user may open the incubation chamber for removal or inclusion of dishes or microtubes at any time. Upon closing the lid, gas exchange and replacement will occur automatically to ensure optimal growing conditions.

Note: Whenever the camera is opened, a message will appear on the display.

CLEANING THE CULTURE CHAMBERS

Use cloth or soft paper towel for cleaning. Gently wipe the incubation chamber with lint free cloth or soft paper towel moistened with 70% alcohol to sanitize.

If a culture medium spill occurs within the culture chamber, remove excess media with a paper towel then clean with a paper towel damp with a solution of water and detergent. Dry completely with paper towels.

Never use abrasive cleaners or solvents as they can damage the metal coating.
Prevent liquids from entering the gas exhaust holes during cleaning.

ALERTS

The equipment will display visual alarms to indicate specific error conditions including: Chamber open; No input pressure; Input pressure out of range; Possible leak; Sensor error; Temperature out of limits; High chamber pressure. In the event of one of these situations, the display will turn red and show the message “ATTENTION”, with the text referring to the problem.


To exit this screen, tap the "X"

USB Reserved for system updates. It is recommended to leave this disabled.

Wi-Fi The "Wi-Fi" mode makes it possible to connect the incubator to the WTA website (<https://wtavet.web.app/br>) and, thus, monitor its performance and adjust parameters including temperature and pressure while monitoring the Datalog from a distance.

Select: "Menu" =>"Operation" =>"Wi-Fi".

To add a WIFI network, enter the wireless network SSID and password. => "Confirm".

The "Wi-Fi" icon will appear in the upper right corner of the display panel. If the "Wi-Fi" device or website does not connect within 1 minute, the mode will be disabled.

To disconnect the equipment from the "Wi-Fi" network:

Select: "Menu" =>"Operation" =>"Wi-Fi" => "Disconnect".

WTA APP FOR REMOTE CONTROL

To use the WTA Website (<https://wtavet.web.app/br>) capabilities, the user must first visit the website to register and accept the Privacy Terms. (Fig. 4-A). An electronic confirmation email will be sent to the selected registration email account. To connect the incubator to the website, enter: “E-mail” and “Password” (Fig 4-B).

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B

Conecte-se

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Figure 4. Registration page (A) and webpage login page (B).

On the home page of the WTA portal (Fig. 5) the user can monitor the environmental parameters of the incubation chambers; access “SETTINGS” to make adjustments or information stored in the “DATALOG”.

EVE18014

	Chamber 1	Chamber 2	Chamber 3
Status	Online	Online	Offline
Plate (°C)	38.5	38.5	38
Lid (°C)	38.6	38.7	38.2
Bubbler (°C)	36	36	38
Setpoint (°C)	38.5	38.5	38
Gas Inlet (Psi)	15	15	0
Chamber (mmHg)	5	5	0
Setpoint (mmHg)	5	5	5
Last Update (HH:mm)	14:31		

SETTINGS

DATALOG

REMOVE

“Online” = Chamber connected to the internet

“Offline” = Chamber not connected to the internet

Figure 5. Home Page.

To make adjustments, press “SETTINGS”. Each chamber can be viewed and programmed individually (Fig. 6). To change the Temperature or Pressure of the chamber, use the slide to select the desired value. Once selections are complete, press the “update” icon

To initiate a gas exchange sequence of the chamber, press the respective update icon

To update the date and time, press the respective update icon




Figure 6. Chamber 1 settings screen.

To access the information stored in the incubator's internal database, press "DATALOG". Then select the current date and press send.

An individual graph will be displayed for each chamber with the Temperature (Fig. 7-A) and Pressure (Fig. 7-B) values for the last 24 hours of operation. To switch between <Temperature> or <Pressure> graphics, press the respective button.

A

B



Figure 7. Chamber 1 datalog graphs.
Temperature (A) and Pressure (B).

The “Datalog Analysis” function shows any error or deviations in chamber conditions (and their time) during the selected time frame.

The Temperature and Pressure data can be exported in PDF format (Fig.8) from the selected date by pressing the “PDF” icon (Fig. 8 -).

EVE18014 - Chamber 1 - 01/06/21

Date	Hour	Plate (°C)	Lid (°C)	Bubbler (°C)	Setpoint (°C)	Chamber (mmHg)	Setpoint (mmHg)	Gas Inlet (Psi)	Status
01/06/21	22:33:43	38.5	38.8	36	38.5	6	5	9	0
01/06/21	22:38:45	38.5	38.8	35.9	38.5	6	5	10	0
01/06/21	22:43:47	38.5	38.7	36	38.5	6	5	10	3
01/06/21	22:48:49	38.5	38.7	36	38.5	6	5	10	9
01/06/21	22:53:51	38.5	38.6	36	38.5	5	5	10	0
01/06/21	22:58:53	38.5	38.7	35.9	38.5	6	5	10	0
01/06/21	23:03:55	38.5	38.8	36	38.5	6	5	10	0
01/06/21	23:13:59	38.5	38.7	35.9	38.5	5	5	10	0
01/06/21	23:19:01	38.5	38.8	36	38.5	6	5	10	0
01/06/21	23:24:03	38.5	38.7	35.9	38.5	6	5	10	0
01/06/21	23:29:05	38.4	38.6	36	38.5	5	5	10	0
01/06/21	23:34:07	38.5	38.6	35.9	38.5	5	5	9	0
01/06/21	23:39:09	38.5	38.7	36	38.5	6	5	8	4
01/06/21	23:44:11	38.4	38.8	35.9	38.5	5	5	8	0
01/06/21	23:49:13	38.5	38.8	36	38.5	6	5	9	0
01/06/21	23:54:15	38.5	38.7	35.9	38.5	6	5	10	0
01/06/21	23:56:59	38.5	38.7	36	38.5	5	5	10	2
01/06/21	23:57:10	38.4	38.8	36	38.5	6	5	4	7
01/06/21	23:59:17	38.5	38.8	36	38.5	8	5	6	0

Figure 8. Chamber 1 Datalog export table (PDF).

This report also shows the gas conditions, any changes designated and errors during operation indicated by a numerical codes displayed in the “Status” column (Fig. 9).

Status	Description
0	Normal
1	Manual gas exchange initiated
2	Automatic gas exchange initiated
3	Closed lid and gas exchange initiated
4	Chamber opened
5	Insufficient inlet pressure
6	No inlet pressure
7	Inlet pressure beyond designated range
8	Possible gas leak
9	Gas exchange completed successfully
10	Pressure adjusted
11	Temperature sensor error
12	Chamber temperature out of range
13	Chamber pressure high
14	Chamber base temperature out of range
15	Chamber lid temperature out of range
55	Chamber turned on

Figure 9. Command and Error Codes

Notifications (“Push Notifications”). After logging in to the WTA web portal, the user may receive error notifications during an incubation session such as: i) Loss/Reestablishment of internet connection (indicating a potential power outage in the laboratory or internet disconnection); ii) Open Chamber; iii) Gas leak; iv) No inlet pressure (indicating either empty gas cylinder or poorly adjusted pressure gauge); v) Pressure outside the set range; vi) Temperature outside the set range.

Note: For this online monitoring to work, it is necessary to log in at least once to your WTA web portal and enable your browser to receive notifications.

WARRANTY

WTA – Watanabe Applied Technology Eireli - EPP follows strict internal quality control standards in the manufacturing process of its products. In addition to the legal warranty of 90 (ninety) days, provided by the Brazilian Consumer Defense Code, WTA guarantees an additional warranty of 33 months from the date of purchase, provided that the equipment has been installed and used in accordance with the guidelines contained within the user's manual.

In cases of ownership transfer, the warranty will continue to be valid according to the terms and conditions specified above.

This product has been carefully evaluated at every step of the manufacturing process. This warranty will always accompany the product, valid only in the country where the product was purchased.

The warranty is understood as the free repair of the product and/or replacement of parts that, according to the technical opinion of the WTA, present a manufacturing defect.

Only WTA Technical Assistance professionals are authorized to repair defects under the warranty, upon presentation of the original sales invoice from the first purchaser.

Included Warranty Services and Coverage

- WTA technical assistance.
- General parts used in maintenance, subject to the exceptions provided for in item "Exclusion of Guarantee".

The product must be sent directly to WTA for repair or return, noting that:

- The product must be properly packaged, if possible in its original packaging, along with the manual, cables and other accessories that have been delivered at the time of purchase.
- The product will be assessed by WTA technicians prior to repair to identify problems.
- Transport, postage and packaging expenses will be the responsibility of WTA during the legal warranty period of 90 days. After this period such costs will be the responsibility of the customer.

Loss of Warranty

This warranty does not cover damage or defects arising from the following situations:

- Installation and use of the product inappropriately or in disagreement with the respective

guidelines and information described in this Manual, including installation, external adjustment, cleaning and use of accessories not in accordance with specifications defined within.

- Defects caused by accidents, falls, natural disaster (rain, flood, lightning strike), excessive humidity (water, liquids in general), salinity, theft, loss or improper transportation of the product.
- Installation, adjustment, repair or technical assistance carried out by people/ parties outside of WTA technical assistance
- Invoice and/or serial number tampered with or removed.
- Modifications, alterations or violations of the original product.
- Product connected to improper electrical network or subject to excessive electrical fluctuations.
- A showcase or demonstration unit is purchased or if the purchase price was lower than market value due to atypical purchase arrangements.

Important: WTA is not responsible for damage caused to equipment during shipping. Upon delivery, ensure and carefully inspect that it is not damaged or defective. If you notice any irregularities, return the product immediately.

The Warranty does not include

- Periodic maintenance
- Travel and/or technical expenses accrued by WTA technical assistance.

Responsibilities

WTA has no control or influence over the conditions under which the equipment is used, or the methodologies employed by users. Thus, WTA does not assume any responsibility for the results, use and/or performance of the equipment (including: damages for lost profits, interruption of commercial activity, loss of commercial information or other monetary losses). WTA expects the equipment to be used by trained personnel. Purchaser acknowledges that the equipment is subject to all terms and conditions outlined in this manual.

COMMON PROBLEMS AND SOLUTIONS

Many errors and functional problems can often be easily fixed. If you notice any problem with the incubator, consider the following potential solutions:

Display is Black, Equipment does Not Turn On

- Check if the equipment is turned on. The “ON/OFF” button will be illuminated when ON.
- Turn off the equipment, wait 5 seconds and turn on again.

Incubator or Incubator Chambers Overheating

- Check that the temperature is set correctly.
- Check if the cables are properly connected.
- Turn off and unplug the incubator, wait 5 seconds and turn on again.

Malfunction of Gas Exchange

- Check that the equipment inlet pressure is at 10 Psi.
- Check if the amount of water in the bubbler is at the recommended level.

No Cellular Connection

- Ensure the incubator is connected to Wi-Fi, the symbol should be visible on the display.
- Disconnect the equipment from the current Wi-Fi network and reconnect it.

“Chamber Open”

- Check that the lid is completely closed.
- Check the condition of the silicone seal on the lid.
- Perform a manual gas exchange to check if the fault persists.

“No Input Pressure”

- Check if the gas hose from the cylinder regulator is properly connected to the equipment.

“Input Pressure Out of Range”

- Adjust the gas inlet regulator until the display reads 10 Psi.

“Possible Leak”

- Check the connection of all hoses.
- Check for debris and seal of the chamber lids.

If the problem persists, contact WTA technical assistance.

Note: it is not possible to connect the equipment to networks that have a proxy.

TECHNICAL ASSISTANCE AND CUSTOMER SERVICE

For questions or additional information about the warranty, assembly, installation, adjustment and use of the product, please contact our customer service department:

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EVE Benchtop Incubator - WTA User Guide.

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Images contained in the Guide are merely illustrative.

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Todos os direitos reservados | All rights reserved